

RAM®

# WHEELCHAIR MOUNT

COMPATIBILITY GUIDE

Phone • Tablet • Camera • Cup Holder



# TABLE OF CONTENTS

<b>2</b>	How it Works
<b>3</b>	Light-Duty vs Heavy-Duty
<b>4</b>	How to Build Your Wheelchair Mount
<b>5</b>	Holders
<b>6</b>	Bases & Arms
<b>7</b>	21st Century Scientific Compatibility Chart
<b>8</b>	Invacare Compatibility Chart
<b>9</b>	Meyra Compatibility Chart
<b>10</b>	Ottobock Compatibility Chart
<b>11</b>	Permobil & Pride Mobility Compatibility Charts
<b>12</b>	Quickie & Zippie Compatibility Charts

# HOW IT WORKS

1. **Base**
2. **Arm**
3. **Holder**



## LIGHT-DUTY VS HEAVY-DUTY



### Light-Duty

#### SMALL DEVICES:

B Size 1" Ball  
Phones, Drinks, Cameras, 7"-8" Tablets



### Heavy-Duty

#### LARGE DEVICES:

C Size 1.5" Ball  
9"-11" Tablets

## HOW TO BUILD YOUR WHEELCHAIR MOUNT

1. Select the proper Holder for the device you are mounting (pg. 5)
2. Select your wheelchair make/model in the Compatibility Chart (pg. 7-12)
3. In the Compatibility Chart, select the proper Base and Arm suitable for your device

### Your Complete

# Kit

will be the Holder from Step 1 paired with the Base and Arm from Step 3


# HOLDERS



**RAM-HOL-UN7BU**  
RAM® X-Grip®  
Phone Holder



**RAM-HOL-UN10BU**  
RAM® X-Grip® Large  
Phone Holder



**RAM-HOL-UN8BU**  
RAM® X-Grip® 7"-8"  
Tablet Holder



**RAM-B-202-UN9U**  
RAM® X-Grip® 9"-11"  
Tablet Holder  
For heavy-duty use only.

## ADDITIONAL NOTES

- All Holders are compatible with all Base/Arm options in this guide
- Large Tablets are recommended only for Heavy-Duty mounting options
- Some chair mounting locations will only recommend Light-Duty or Heavy-Duty
- The Leg Rail Mounting location is suitable for both Light-Duty and Heavy-Duty use



**RAP-B-202U-GOP1**  
Action Camera  
Ball Adapter



**RAP-B-417BU**  
RAM® Level-Cup™ XL  
32oz Cup Holder

# BASES & ARMS

## LIGHT-DUTY



**RAM-B-238-WCT**  
Wheelchair  
Track Base



**RAM-B-400U**  
RAM® Tough-Claw™  
Rail Base



**RAM-B-201U-C**  
6" Socket Arm for  
Arm Rests



**RAM-B-201-201U-C**  
10" Socket Arm for  
Arm Rests



**RAM-B-200-9-201**  
16" Pipe & Socket Arm  
for Seat Tracks

## HEAVY-DUTY



**RAP-400-18-B-201**  
RAM® Tough-Claw™  
with 26" Extension Arm



**RAM-200-9-BC-201**  
19" Pipe & Socket  
Arm for Seat Tracks



**RAM-238-WCT**  
Wheelchair  
Track Base

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
Boulder Series	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Seat Rail	Base and Arm: RAP-400-18-B-201
Boulder 300M	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Seat Rail	Base and Arm: RAP-400-18-B-201
Boulder Plus HVLP	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Seat Rail	Base and Arm: RAP-400-18-B-201
Boulder Plus Narrow Base	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Seat Rail	Base and Arm: RAP-400-18-B-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6	ADDITIONAL NOTES
Storm Series Torque 3	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
Storm Series Torque 3 (Power Seating)	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
	Seat Rail	Base and Arm: RAP-400-18-B-201	
Storm Series Torque 3G (Power Seating)	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
TDX SP (Power Seating)	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
TDX SP2	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
TDX SP2 (Single Power Seating)	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
TDX SP2 (Multiple Power Seating)	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C	Base connects with provided metric screws & nuts
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
TDX SP2X (Single Power Seating)	Arm Rail	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
TDX SP2X (Multiple Power Seating)	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C	Base connects with provided metric screws & nuts
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
Avanti 1.736	Leg Rail	Base and Arm: RAP-400-18-B-201
Avanti Pro 1.735	Leg Rail	Base and Arm: RAP-400-18-B-201
Brix 1.123	Leg Rail	Base and Arm: RAP-400-18-B-201
Budget 9.050	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair 1.750	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair 1.850	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair 2 Polaro 2.845	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair2 HD 2.865	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair2 Pro 2.850	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair2 Stock Version 2.750	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair2 XXL 2.850	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair Hemi-Special 1.840	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair Lagervariante 1.750	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Eurochair XXL 1.760	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
Flash 1.135	Leg Rail	Base and Arm: RAP-400-18-B-201
Format 3.940	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
HCDA mobile commode XXL	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Seat Rail	Base and Arm: RAP-400-18-B-201
	Leg Rail	Base and Arm: RAP-400-18-B-201
Hurricane Assist	Leg Rail	Base and Arm: RAP-400-18-B-201
Hurricane Everyday 1.880	Leg Rail	Base and Arm: RAP-400-18-B-201
Hurricane Pro	Leg Rail	Base and Arm: RAP-400-18-B-201
Hurricane S 1.880	Leg Rail	Base and Arm: RAP-400-18-B-201
Hurricane Sport 1.880	Leg Rail	Base and Arm: RAP-400-18-B-201
Motivo 2.250	Leg Rail	Base and Arm: RAP-400-18-B-201
Service 3.600	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
Service 3.600	Leg Rail	Base and Arm: RAP-400-18-B-201
SMART F 2.360	Leg Rail	Base and Arm: RAP-400-18-B-201
SMART S 2.370	Leg Rail	Base and Arm: RAP-400-18-B-201
Solero Light 9.072	Leg Rail	Base and Arm: RAP-400-18-B-201
XR 1.911	Leg Rail	Base and Arm: RAP-400-18-B-201
XStar 1.160	Leg Rail	Base and Arm: RAP-400-18-B-201
ZX1 1.360	Leg Rail	Base and Arm: RAP-400-18-B-201
ZX3 1.370	Leg Rail	Base and Arm: RAP-400-18-B-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
Avantgarde DS	Leg Rail	Base and Arm: RAP-400-18-B-201
Avantgarde DV	Leg Rail	Base and Arm: RAP-400-18-B-201
Motus	Leg Rail	Base and Arm: RAP-400-18-B-201
Ventus	Leg Rail	Base and Arm: RAP-400-18-B-201
Voyager Evo	Leg Rail	Base and Arm: RAP-400-18-B-201
Zenit	Leg Rail	Base and Arm: RAP-400-18-B-201
Pointer	Leg Rail	Base and Arm: RAP-400-18-B-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
F3 Corpus	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
F5 Corpus	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
F5 Corpus VS	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
K300 PS Jr.	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
K450 MX	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
Jazzy Passport	Seat Rail	Base and Arm: RAP-400-18-B-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
M1	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
M3 Corpus	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
M300 Corpus HD	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
M300 PS Jr.	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
M5 Corpus	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
TiLite Series	Leg Rail	Base and Arm: RAP-400-18-B-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
5R+B102:E135	Leg Rail	Base and Arm: RAP-400-18-B-201
7 Series	Leg Rail	Base and Arm: RAP-400-18-B-201
All Court	Leg Rail	Base and Arm: RAP-400-18-B-201
GP Series	Leg Rail	Base and Arm: RAP-400-18-B-201
IRIS	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
M6	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Match Point	Leg Rail	Base and Arm: RAP-400-18-B-201
P-222 SE	Arm Rail	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Pulse Series	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201
Q7	Leg Rail	Base and Arm: RAP-400-18-B-201
Q700 M	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
GS	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
IRIS	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
Kidz	Leg Rail	Base and Arm: RAP-400-18-B-201

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6	ADDITIONAL NOTES
QM-7	Arm Rest	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-201-201U-C	Base connects with provided metric screws & nuts
QM-7	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
QRi	Leg Rail	Base and Arm: RAP-400-18-B-201	
Quickie 2 Family	Leg Rail	Base and Arm: RAP-400-18-B-201	
QXi/QX	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Leg Rail	Base and Arm: RAP-400-18-B-201	
S-6 Series	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Seat Track	Light-Duty Base: RAM-B-238-WCT Light-Duty Arm: RAM-B-200-9-201 Heavy-Duty Base: RAM-238-WCT Heavy-Duty Arm: RAM-200-9-BC-201	
		Leg Rail	Base and Arm: RAP-400-18-B-201
SR45	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Leg Rail	Base and Arm: RAP-400-18-B-201	
Xenon2 Series	Leg Rail	Base and Arm: RAP-400-18-B-201	
Xtender	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C	
	Seat Rail	Base and Arm: RAP-400-18-B-201	
	Leg Rail	Base and Arm: RAP-400-18-B-201	

CHAIR MODEL	MOUNTING LOCATIONS	RAM® PART NUMBER(S) See pg. 6
TS	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201
X-CAPE	Leg Rail	Base and Arm: RAP-400-18-B-201
Zippie 2	Arm Rest	Light-Duty Base: RAP-B-400U Light-Duty Arm: RAM-B-201U-C
	Leg Rail	Base and Arm: RAP-400-18-B-201